

Anti-Lipocalin-2/NGAL (3C6) Mouse antibody



Catalog Number: 168265

Product name

Monoclonal

Anti-Lipocalin-2/NGAL (3C6) Mouse antibody

Ig Type

Mouse IgG2a

Specificity

Human

Applications

Antibody description

ELISA

Lipocalin-2/NGAL (3C6) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human NGAL protein expressed in E.coli.

Formulation


Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300.

Storage

Store at 4°C short term.Avoid freeze / thaw cycle.

Clonality

Validations



Anti-FITC mouse mAb was used to capture FITC-labeled NGAL Mouse mAb (168265), HRP-labeled NGAL Mouse mAb (201179-2F6) was used for detecting.