

CD46 Polyclonal Antibody



Catalog Number: 161665

Product name

CD46 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to CD46

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 40-320 of human CD46 (NP_722548.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

Validations



Western blot - CD46 Polyclonal Antibody

Western blot analysis of extracts of HT-29 cells, using CD46 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunofluorescence - CD46 Polyclonal Antibody

Immunofluorescence analysis of HepG2 cells using CD46 antibody . Blue: DAPI for nuclear staining.