

# Recombinant RET4 protein (His tag)



Catalog Number: 546665

---

## Gene Name Synonym

Plasma retinol-binding protein; PRBP; RBP; 1-182; 1-181; 1-179; 1-176

## Protein Construction

DNA encoding 19-201 aa (P02753) were fused with 6His tag.

## Organism

Human

## Expression Host

E. coli

## Activity

Not tested.

## Purity

>87%, by SDS-PAGE with Coomassie Brilliant Blue staining.

## Endotoxin

Please contact the lab for more information.

## Stability

Samples are stable for up to twelve months from

date of receipt at -70°C

## Molecular Mass

Predictes a molecular mass of 22.73 kDa. In SDS-PAGE under reducing conditions, it migrates as an approximately 24 kDa band.

## Formulation

Lyophilized from sterile PBS, pH 7.4

## Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

## Reconstitution

Reconstitute at 0.25 µg/µl in sterile water for short-term storage. Reconstitution with 50% glycerol solution is recommended for longer term storage (see Stability and Storage for more details). If a different concentration is needed for your purposes please adjust the reconstitution volume as required (please note: the ion concentration of the final solution will vary according to the volume used). Note: Centrifuge vial before opening. When reconstituting, gently pipet and wash down the sides of the vial to ensure full recovery of the protein into solution.