

# TSPAN7 antibody



Catalog Number: 116365

## Product name

TSPAN7 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

TSPAN7 Rabbit Polyclonal antibody. Positive IHC detected in human spleen tissue. Positive IP detected in K-562 cells. Positive WB detected in K-562 cells. Observed molecular weight by Western-blot: 28-35kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_004615). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC, IP

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000


IP: 1:200-1:2000

IHC: 1:20-1:200


## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:116365(TSPAN7 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human spleen using Catalog No:116365(TSPAN7 antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemical of paraffin-embedded human spleen using Catalog No:116365(TSPAN7 antibody) at dilution of 1:50 (under 40x lens)

# TSPAN7 antibody

Catalog Number: 116365

---



IP Result of anti-TSPAN7 (IP:Catalog No:116365, 4ug; Detection:Catalog No:116365 1:500) with K-562 cells lysate 3600ug.