

BAX antibody



Catalog Number: 108365

Product name

BAX antibody

WB: 1:500-1:5000

IP: 1:500-1:5000

Specificity

Human, Mouse, Monkey; other species not tested.

IHC: 1:20-1:200

Antibody description

BAX Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, COLO 320 cells, COS-7 cells, HeLa cells, MCF7 cells, RAW264.7 cells. Positive IP detected in HEK-293 cells. Positive IHC detected in human kidney tissue, human breast cancer tissue. Observed molecular weight by Western-blot: 21 kDa

Preparation

This antibody was obtained by immunization of BAX recombinant protein (Accession Number: NM_138761). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IP


Dilutions

Recommended Dilution:


Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108365(BAX antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:108365(BAX Antibody) at dilution of 1:50 (under 10x lens)




BAX antibody



Catalog Number: 108365

Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:108365(BAX Antibody) at dilution of 1:50

IP Result of anti-BAX (IP:Catalog No:108365, 4ug; Detection:Catalog No:108365 1:1000) with HEK-293 cells lysate 1480ug.



(under 40x lens)