

PP11 antibody



Catalog Number: 114075

Product name

PP11 antibody

IHC: 1:20-1:200

IF: 1:50-1:500

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PP11 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive IHC detected in human breast cancer tissue. Positive WB detected in mouse lung tissue, HeLa cells. Observed molecular weight by Western-blot: 35-40 kDa

Preparation

This antibody was obtained by immunization of PP11 recombinant protein (Accession Number: NM_006025). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:114075(P11 antibody) at dilution of 1:600



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:114075(P11 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded

PP11 antibody



Catalog Number: 114075

human breast cancer using Catalog

No:114075(P11 antibody) at dilution of 1:100
(under 40x lens)