

ANKRD33 antibody



Catalog Number: 108075

Product name

ANKRD33 antibody

Specificity

Human; other species not tested.

Antibody description

ANKRD33 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue. Positive WB detected in Y79 cells. Observed molecular weight by Western-blot: 49 kDa

Preparation

This antibody was obtained by immunization of ANKRD33 recombinant protein (Accession Number: NM_001130015). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



Y79 cells were subjected to SDS PAGE followed by western blot with Catalog No:108075(ANKRD33 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:108075(ANKRD33 Antibody) at dilution of 1:50 (under 40x lens)