

Gastrin antibody

Catalog Number: 110875

Product name

Gastrin antibody

Specificity

Human; other species not tested.

Antibody description

Gastrin Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells. Positive IHC detected in human stomach tissue.

Preparation

This antibody was obtained by immunization of Gastrin recombinant protein (Accession Number: NM_000805). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human stomach using Catalog No:110875(GAST antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human stomach using Catalog No:110875(GAST antibody) at dilution of 1:50 (under 40x lens)



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110875(GAST Antibody) at dilution of 1:500