

# COG2 Polyclonal Antibody



Catalog Number: 161175

---

## Product name

COG2 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to COG2

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human COG2 (NP\_031383.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IF

## Dilutions

IF 1:20 - 1:100

## Validations



Immunofluorescence - COG2 Polyclonal Antibody

Immunofluorescence analysis of MCF-7 cells using COG2 antibody . Blue: DAPI for nuclear staining.