

Anti-CDKN2A/p16INK4a Rabbit antibody



Catalog Number: 167175

Product name

Anti-CDKN2A/p16INK4a Rabbit antibody

Specificity

Human

Antibody description

CDKN2A/p16INK4a Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human CDKN2A/p16INK4a

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC


Dilutions

WB: 1/1000

ICC/IF: 1/20-1/100

FC: 1/20



Validations



Immunocytochemistry of CDKN2A/p16INK4a(green) in HeLa cells using CDKN2A/p16INK4a Rabbit pAb at dilution 1/50,

and DAPI(blue)

Western blot detection of CDKN2A/p16INK4a in 293 cell lysates using CDKN2A/p16INK4a Rabbit pAb(1:1000 diluted).Predicted band size:17KDa.Observed band size:17KDa.



Western blot detection of CDKN2A-p16INK4a in 3T3,HeLa cell lysates using CDKN2A/p16INK4a Rabbit pAb(1:1000 diluted).Predicted band size:17KDa.Observed band size:17KDa.