

Product name

CRNN antibody

IHC: 1:20-1:200

Specificity

Human, Mouse; other species not tested.

Antibody description

CRNN Rabbit Polyclonal antibody. Positive IHC detected in human oesophagus tissue, human cervical cancer tissue, human skin tissue. Positive WB detected in mouse brain tissue, HepG2 cells. Observed molecular weight by Western-blot: 70 kDa, 58 kDa

Preparation

This antibody was obtained by immunization of CRNN recombinant protein (Accession Number: NM_016190). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

Validations

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109575(CRNN antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human oesophagus tissue slide using Catalog No:109575(CRNN Antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieved with Citric acid buffer, pH6.0



Immunohistochemistry of paraffin-embedded

human oesophagus tissue slide using Catalog
No:109575(CRNN Antibody) at dilution of 1:200

(under 40x lens). Heat mediated antigen retrieved
with Citric acid buffer, pH6.0