

TADA2L antibody



Catalog Number: 115975

Product name

TADA2L antibody

IP: 1:200-1:1000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TADA2L Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Positive IP detected in mouse testis tissue. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 51.5 kDa

Preparation

This antibody was obtained by immunization of TADA2L recombinant protein (Accession Number: NM_001488). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:115975(TADA2L antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:115975(TADA2L antibody) at dilution of 1:50 (under 10x lens)



IP Result of anti-TADA2L (IP:Catalog No:115975,

TADA2L antibody



Catalog Number: 115975

4ug; Detection:Catalog No:115975 1:300) with mouse testis tissue lysate 4000ug.