

MLF2 antibody



Catalog Number: 112675

Product name

MLF2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MLF2 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Positive IHC detected in human gliomas tissue. Observed molecular weight by Western-blot: 28 kDa

Preparation

This antibody was obtained by immunization of MLF2 recombinant protein (Accession Number: NM_001382226). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112675(MLF2 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:112675(MLF2 antibody) at dilution of 1:50 (under 10x lens)