

SF3B2 Polyclonal Antibody



Catalog Number: 164675

Product name

SF3B2 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to SF3B2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 616-895 of human SF3B2 (NP_006833.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations




Western blot - SF3B2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using SF3B2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.



Immunohistochemistry - SF3B2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat liver using SF3B2 antibody at dilution of 1:100




(40x lens).

SF3B2 Polyclonal Antibody

Catalog Number: 164675

Immunohistochemistry - SF3B2 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat kidney using SF3B2 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - SF3B2 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat heart using SF3B2 antibody at dilution of 1:100



(40x lens).


Immunohistochemistry - SF3B2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human lung cancer using SF3B2 antibody at dilution of 1:100 (40x lens).




Immunohistochemistry - SF3B2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver cancer using SF3B2 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - SF3B2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human oophoroma using SF3B2 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - SF3B2 Polyclonal Antibody


SF3B2 Polyclonal Antibody



Catalog Number: 164675


Antibody

Immunohistochemistry of paraffin-embedded human kidney cancer using SF3B2 antibody at dilution of 1:100 (40x lens).




Immunohistochemistry - SF3B2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human gastric cancer using SF3B2 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - SF3B2 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded mouse Intestine using SF3B2 antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - SF3B2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse liver using SF3B2 antibody at dilution of 1:100 (40x lens).



Immunofluorescence - SF3B2 Polyclonal Antibody

Immunofluorescence analysis of A-549 cells using SF3B2 antibody .