

PSMA4 antibody



Catalog Number: 114375

Product name

PSMA4 antibody

Immunogen

[Human PSMA4 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PSMA4 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, BxPC-3 cells, mouse kidney tissue, mouse spleen tissue, U-937 cells. Positive IP detected in BxPC-3 cells. Observed molecular weight by Western-blot: 29 kDa

Preparation

This antibody was obtained by immunization of PSMA4 recombinant protein (Accession Number: BC056249). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:114375(PSMA4 Antibody) at dilution of 1:1000



IP Result of anti-PSMA4 (IP:Catalog No:114375, 4ug; Detection:Catalog No:114375 1:500) with BxPC-3 cells lysate 2800ug.