

TTBK2 antibody



Catalog Number: 116375

Product name

TTBK2 antibody

IHC: 1:20-1:200

IF: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TTBK2 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, SH-SY5Y cells. Positive IHC detected in human brain tissue, human prostate cancer tissue. Positive IF detected in SH-SY5Y cells. Observed molecular weight by Western-blot: 137-150 kDa

Preparation

This antibody was obtained by immunization of TTBK2 recombinant protein (Accession Number: NM_173500). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:116375(TTBK2 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human brain using Catalog No:116375(TTBK2 antibody) at dilution of 1:100 (under 10x lens)




Immunohistochemical of paraffin-embedded human brain using Catalog No:116375(TTBK2 antibody) at dilution of 1:100 (under 10x lens)

TTBK2 antibody



Catalog Number: 116375

antibody) at dilution of 1:100 (under 40x lens)



Immunofluorescent analysis of (-20oc Ethanol)
fixed SH-SY5Y cells using Catalog
No:116375(TTBK2 Antibody) at dilution of 1:50
and Alexa Fluor 488-conjugated AffiniPure Goat
Anti-Rabbit IgG(H+L)