

# ANKS1B antibody



Catalog Number: 108085

## Product name

ANKS1B antibody

IP: 1:500-1:5000

IHC: 1:20-1:200

## Specificity

Human, Mouse; other species not tested.

IF: 1:20-1:200

## Antibody description

ANKS1B Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in A549 cells, HeLa cells, mouse testis tissue. Positive IHC detected in human testis tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 65-70 kDa

## Preparation

This antibody was obtained by immunization of ANKS1B recombinant protein (Accession Number: NM\_181670). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IF, IHC, IP

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000


## Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:108085(ANKS1B Antibody) at dilution of 1:1000



IP Result of anti-ANKS1B (IP:Catalog No:108085, 4ug; Detection:Catalog No:108085 1:1000) with HeLa cells lysate 3200ug.




# ANKS1B antibody

Catalog Number: 108085


---

Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:108085(ANKS1B Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded

human testis tissue slide using Catalog No:108085(ANKS1B Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells using Catalog No:108085(ANKS1B Antibody) at dilution of 1:50 and Rhodamine-Goat anti-Rabbit IgG