

ARHGEF18 antibody



Catalog Number: 108185

Product name

ARHGEF18 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ARHGEF18 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Positive IHC detected in human pancreas cancer tissue, human kidney tissue. Observed molecular weight by Western-blot: 120-130 kDa

Preparation

This antibody was obtained by immunization of ARHGEF18 recombinant protein (Accession Number: BC014994). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:108185(ARHGEF18 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:108185(ARHGEF18 antibody) at dilution of 1:50 (under 10x lens)