

Parkin antibody



Catalog Number: 113585

Product name

Parkin antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Parkin Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, C6 cells, SH-SY5Y cells. Positive IF detected in SH-SY5Y cells. Observed molecular weight by Western-blot: 52 kDa

Preparation

This antibody was obtained by immunization of Parkin recombinant protein (Accession Number: BC022014). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:20-1:200

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113585(Parkin Antibody) at dilution of 1:300



Immunofluorescent analysis of SH-SY5Y cells using Catalog No:113585(Parkin Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)