

Anti-IRAK4 Rabbit antibody



Catalog Number: 168985

Product name

Anti-IRAK4 Rabbit antibody

Specificity

Human

Antibody description

IRAK4 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human IRAK4

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC

Dilutions


WB: 1/1000-1/5000

IHC: 1/20-1/500

ICC/IF: 1/20

FC: 1/20

Validations



Western blot detection of IRAK4 in Jurkat, HeLa cell lysates using IRAK4 Rabbit mAb(1:1000 diluted). Predicted band size: 52kDa. Observed band size: 52kDa.