

Product name

ELN antibody

Specificity

Human; other species not tested.

Antibody description

ELN Rabbit Polyclonal antibody. Positive IHC detected in mouse lung tissue.

Preparation

This antibody was obtained by immunization of ELN recombinant protein (Accession Number: NM_001278916). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC


Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations

Immunohistochemistry of paraffin-embedded mouse lung tissue slide using Catalog No:110285(ELN Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemistry of paraffin-embedded mouse lung tissue slide using Catalog No:110285(ELN Antibody) at dilution of 1:200 (under 40x lens).