

HCN2 antibody



Catalog Number: 111285

Product name

HCN2 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

HCN2 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue, Jurkat cells. Positive IHC detected in human heart tissue, human brain tissue. Observed molecular weight by Western-blot: 97-110 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001194). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IP, IHC

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111285(HCN2 antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human heart using Catalog No:111285(HCN2 antibody) at dilution of 1:50 (under 40x lens)



IP Result of anti-HCN2 (IP:Catalog No:111285, 4ug; Detection:Catalog No:111285 1:800) with

HCN2 antibody



Catalog Number: 111285

mouse brain tissue lysate 4000ug.