

Anti-RAB7 Rabbit antibody



Catalog Number: 167285

Product name

Anti-RAB7 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

RAB7 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human RAB7

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC

Dilutions


WB: 1/1000-1/5000

IHC: 1/20-1/50


ICC/IF: 1/20-1/100

FC: 1/20


Validations



Immunohistochemistry of RAB7 in paraffin-embedded Human lung cancer tissue using RAB7 Rabbit pAb at dilution 1/20



Immunocytochemistry of RAB7(green) in HeLa cells using RAB7 Rabbit pAb at dilution 1/50, and



DAPI(blue)

Western blot detection of RAB7 in K562,C6,3T3,HeLa cell lysates using RAB7 Rabbit pAb(1:1000 diluted).Predicted band size:24kDa.Observed band size:24kDa.