

TAC1 Polyclonal Antibody



Catalog Number: 161685

Product name

TAC1 Polyclonal Antibody

Ig Type

Rabbit IgG

Specificity

Human, Rat

Applications

IF

Antibody description

Polyclonal antibody to TAC1


Dilutions

IF 1:50 - 1:200

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-129 of human TAC1 (NP_003173.1).

Validations



Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Immunofluorescence - TAC1 Polyclonal Antibody

Clonality

Polyclonal

Immunofluorescence analysis of C6 cells using TAC1 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.