

ANAPC2 Polyclonal Antibody



Catalog Number: 165685

Product name

Rabbit IgG

ANAPC2 Polyclonal Antibody

Applications

Specificity

WB

Human, Mouse, Rat

Dilutions

Antibody description

WB 1:500 - 1:2000

Polyclonal antibody to ANAPC2

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 400-500 of human ANAPC2 (NP_037498.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type



Western blot - ANAPC2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ANAPC2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.