

# Anti-CD63 Rabbit antibody



Catalog Number: 166785

## Product name

Rabbit IgG

Anti-CD63 Rabbit antibody

## Applications

## Specificity

WB, IHC-P, ICC/IF

Mouse

## Dilutions

## Antibody description

WB: 1/1000

CD63 Rabbit polyclonal antibody

IHC: 1/20

## Preparation

ICC/IF: 1/20


Antigen: A synthetic peptide of mouse CD63

## Validations

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

RAW264.7



## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

Western blot detection of CD63 in Raw264.7 cell lysates using CD63 Rabbit pAb(1:1000 diluted).Predicted band size:26KDa.Observed band size:40-70KDa.

## Ig Type