

CDC23 antibody



Catalog Number: 109095

Product name

CDC23 antibody

Immunogen

[Human CDC23 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CDC23 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, HeLa cells. Observed molecular weight by Western-blot: 49 kDa

Preparation

This antibody was obtained by immunization of CDC23 recombinant protein (Accession Number: BC005258). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:109095(CDC23 Antibody) at dilution of 1:600