

GBP6 antibody



Catalog Number: 110895

Product name

GBP6 antibody

Immunogen

[Human GBP6 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

GBP6 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Positive IHC detected in human testis tissue, human lung cancer tissue. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of GBP6 recombinant protein (Accession Number: XM_011540835). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110895(GBP6 Antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:110895(GBP6 Antibody) at dilution of 1:50 (under 10x lens)



GBP6 antibody



Catalog Number: 110895

Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog

No:110895(GBP6 Antibody) at dilution of 1:50 (under 40x lens)