

RPL39 antibody



Catalog Number: 114895

Product name

RPL39 antibody

IHC: 1:20-1:200

IF: 1:50-1:500

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RPL39 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, human liver tissue. Positive IHC detected in human liver cancer tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 6 kDa

Preparation

This antibody was obtained by immunization of RPL39 recombinant protein (Accession Number: NM_001000). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:114895(RPL39 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:114895(RPL39 antibody) at dilution of 1:100 (under 10x lens)




Immunohistochemical of paraffin-embedded

RPL39 antibody



Catalog Number: 114895

human liver cancer using Catalog
No:114895(RPL39 antibody) at dilution of 1:100
(under 40x lens)



Immunofluorescent analysis of (10%
Formaldehyde) fixed HepG2 cells using Catalog
No:114895(RPL39 Antibody) at dilution of 1:50
and Alexa Fluor 488-conjugated AffiniPure Goat
Anti-Rabbit IgG(H+L)