

Product name

Polyclonal

ICOS antibody

Ig Type**Specificity**

Rabbit IgG

Human; other species not tested.

Applications**Antibody description**

ELISA, WB

ICOS Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HL-60 cells. Observed molecular weight by Western-blot: 28 kDa, 45-50 kDa


Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Preparation

This antibody was obtained by immunization of ICOS recombinant protein (Accession Number: NM_012092). Purification method: Antigen affinity purified.

Validations**Formulation**

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:111595(ICOS Antibody) at dilution of 1:300

Clonality