

XPO7 antibody



Catalog Number: 116995

Product name

XPO7 antibody

Immunogen

[Human XPO7 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

XPO7 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, human brain tissue, human testis tissue. Observed molecular weight by Western-blot: 100-110kd

Preparation

This antibody was obtained by immunization of XPO7 recombinant protein (Accession Number: NM_015024). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:116995(XPO7 antibody) at dilution of 1:300