

Aminoacylase 1 antibody



Catalog Number: 107995

Product name

Aminoacylase 1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Aminoacylase 1 Rabbit Polyclonal antibody.
Positive IHC detected in human prostate cancer tissue. Positive WB detected in K-562 cells.
Observed molecular weight by Western-blot: 46 kDa

Preparation

This antibody was obtained by immunization of Aminoacylase 1 recombinant protein (Accession Number: NM_000666). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:107995(ACY1 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human prostate cancer tissue slide using Catalog No:107995(ACY1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human prostate cancer tissue slide using Catalog No:107995(ACY1 Antibody) at dilution of 1:50

Aminoacylase 1 antibody



Catalog Number: 107995

(under 40x lens)