

Anti-LYVE1 Rabbit antibody



Catalog Number: 166795

Product name

Rabbit IgG

Anti-LYVE1 Rabbit antibody

Applications

Specificity

WB, IHC-P, FC, IP

Mouse

Dilutions

Antibody description

WB: 1/1000

LYVE1 Rabbit polyclonal antibody

IHC: 1/200

Preparation

FC: 1/20

Antigen: Recombinant protein of mouse LYVE1

IP: 1/20

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage


Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Validations



Western blot detection of LYVE1 in Mouse Lung lysates using LYVE1 Rabbit pAb(1:1000 diluted). Predicted band size: 35KDa. Observed band size: 35-70KDa.