

TFDP1 antibody



Catalog Number: 116006

Product name

TFDP1 antibody

Immunogen

[Human TFDP1 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TFDP1 Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue.

Preparation

This antibody was obtained by immunization of TFDP1 recombinant protein (Accession Number: XM_047430563). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC

Dilutions

Recommended Dilution:
IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:116006(TFDP1 antibody) at dilution of 1:50 (under 10x lens)