

Finotonlimab (Synonyms: SCT I10A [WHO-DD])



Catalog Number: 177006

Product name

Finotonlimab (Synonyms: SCT I10A [WHO-DD])

Specificity

Human PDCD1 / PD-1 / CD279

Antibody description

Anti-PDCD1 / PD-1 / CD279 Antibody
(Finotonlimab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG4SP

Applications


ELISA

Validations



SDS-PAGE

Anti-PDCD1 / PD-1 / CD279 Antibody (finotonlimab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



SEC-HPLC

The purity of Anti-PDCD1 / PD-1 / CD279 Antibody (finotonlimab) is more than 95%, determined by SEC-HPLC.