

Anti-CREB Rabbit antibody



Catalog Number: 166806

Product name

WB: 1/1000

Anti-CREB Rabbit antibody

IHC: 1/20

Specificity

FC: 1/20

Human, Mouse, Rat

IP: 1/20

Antibody description

CREB Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human CREB

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


WB, IHC-P, FC, IP

Dilutions

Validations



Immunohistochemistry of CREB in paraffin-embedded Human lung cancer tissue using CREB



Rabbit pAb at dilution 1/20

Western blot detection of CREB in HeLa cell lysates using CREB Rabbit pAb(1:1000 diluted). Predicted band size:37KDa.Observed band size:43KDa.