

C9orf9 antibody



Catalog Number: 108806

Product name

C9orf9 antibody

Immunogen

[Human C9orf9 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

C9orf9 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, mouse liver tissue, rat testis tissue. Positive IF detected in A431 cells. Observed molecular weight by Western-blot: 18, 30 kDa

Preparation

This antibody was obtained by immunization of C9orf9 recombinant protein (Accession Number: NM_001316897). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:10-1:100

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:108806(C9orf9 Antibody) at dilution of 1:300



Immunofluorescent analysis of A431 cells using Catalog No:108806(C9orf9 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)