

ANT2 antibody



Catalog Number: 108106

Product name

ANT2 antibody

WB: 1:200-1:1000

IHC: 1:20-1:200

Immunogen

[Human SLC25A5 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ANT2 Rabbit Polyclonal antibody. Positive IHC detected in human kidney tissue. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 30-32 kDa

Preparation

This antibody was obtained by immunization of ANT2 recombinant protein (Accession Number: NM_001152). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


Validations



Immunohistochemistry of paraffin-embedded human kidney slide using Catalog No:108106(SLC25A5 Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human kidney slide using Catalog No:108106(SLC25A5 Antibody) at dilution of 1:50



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog

ANT2 antibody



Catalog Number: 108106

No:108106(SLC25A5 Antibody) at dilution of 1:300