

EXOSC1 antibody



Catalog Number: 110506

Product name

EXOSC1 antibody

Immunogen

[Human EXOSC1 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

EXOSC1 Rabbit Polyclonal antibody. Positive WB detected in SGC-7901 cells, HeLa cells, human stomach tissue, Jurkat cells. Positive IHC detected in human lymphoma tissue. Observed molecular weight by Western-blot: 21 kDa

Preparation

This antibody was obtained by immunization of EXOSC1 recombinant protein (Accession Number: NM_016046). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



SGC-7901 cells were subjected to SDS PAGE followed by western blot with Catalog No:110506(EXOSC1 antibody) at dilution of 1:200



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:110506(EXOSC1 antibody) at dilution of 1:100 (under 10x lens)