

Anti-Rad51 Rabbit antibody



Catalog Number: 167906

Product name

Anti-Rad51 Rabbit antibody

Specificity

Human

Antibody description

Rad51 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Rad51

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

IHC-P, ICC/IF, FC


Dilutions

IHC: 1/50

ICC/IF: 1/20

FC: 1/50

Validations



Western blot detection of Rad51 in A549,HL-60,U2OS cell lysates using Rad51 Rabbit mAb(1:1000 diluted).Predicted band size:37kDa.Observed band size:37kDa.