

FcRn-Specific antibody



Catalog Number: 110606

Product name

Polyclonal

FcRn-Specific antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

FcRn-Specific Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue, human liver tissue. Observed molecular weight by Western-blot: 35 kDa

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_004107).
Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:110606(FcGRT-Specific antibody) at dilution of 1:300