

C16orf61 antibody



Catalog Number: 108606

Product name

C16orf61 antibody

Clonality

Polyclonal

Immunogen

[Human C16orf61 Recombinant protein \(GST tag & His tag\)](#)

Ig Type

Rabbit IgG

Specificity

Human; other species not tested.

Applications

ELISA, WB

Antibody description

C16orf61 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, human kidney tissue. Observed molecular weight by Western-blot: 9 kDa

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of C16orf61 recombinant protein (Accession Number: NM_020188). Purification method: Antigen affinity purified.


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:108606(C16orf61 Antibody) at dilution of 1:600