

# ZNF649 antibody



Catalog Number: 117016

## Product name

ZNF649 antibody

## Immunogen

[Human ZNF649 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ZNF649 Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells, mouse heart tissue. Observed molecular weight by Western-blot: 50-57 kDa

## Preparation

This antibody was obtained by immunization of ZNF649 recombinant protein (Accession Number: NM\_023074). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:117016(ZNF649 Antibody) at dilution of 1:300