

Anti-CACNA1H antibody



Catalog Number: 175816

Product name

Anti-CACNA1H antibody

Specificity

Human, Rat, Mouse

Antibody description

Rabbit polyclonal antibody to CACNA1H

Preparation

This antigen of this antibody was synthetic peptide of cav3.2

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

IHC-P

Dilutions

IHC-P: 1:50-100

Validations




Fig1.; Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue using Cav3.2Rabbit pAb diluted at 1:200.