

SUB1 antibody



Catalog Number: 115816

Product name

SUB1 antibody

IHC: 1:20-1:200

IF: 1:50-1:500

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SUB1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, A549 cells. Positive IHC detected in human pancreas cancer tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 20 kDa, 35 kDa

Preparation

This antibody was obtained by immunization of SUB1 recombinant protein (Accession Number: XM_047416661). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115816(SUB1 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:115816(SUB1 antibody) at dilution of 1:50 (under 25x lens)



Immunofluorescent analysis of (10%

SUB1 antibody



Catalog Number: 115816

Formaldehyde) fixed HeLa cells using Catalog No:115816(SUB1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)