

# CYTSB antibody



Catalog Number: 109816

## Product name

CYTSB antibody

## Immunogen

[Human CYTSB Recombinant protein \(GST tag & His tag\)](#)

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

CYTSB Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, HeLa cells. Observed molecular weight by Western-blot: 119kd

## Preparation

This antibody was obtained by immunization of CYTSB recombinant protein (Accession Number: NM\_001243438). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:109816(CYTSB antibody) at dilution of 1:1000