

KPNA6 antibody



Catalog Number: 112116

Product name

KPNA6 antibody

Immunogen

[Human KPNA6 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

KPNA6 Rabbit Polyclonal antibody. Positive WB detected in BxPC-3 cells, HeLa cells. Positive IHC detected in human prostate cancer tissue. Observed molecular weight by Western-blot: 55-60 kDa

Preparation

This antibody was obtained by immunization of KPNA6 recombinant protein (Accession Number: NM_012316). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:112116(KPNA6 antibody) at dilution of 1:50



(under 10x lens)

BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:112116(KPNA6 Antibody) at dilution of 1:600