

Anti-FEN1 (7H8) Mouse antibody



Catalog Number: 168116

Product name

Anti-FEN1 (7H8) Mouse antibody

Specificity

Human, Mouse, Rat, Monkey, Hamster

Antibody description

FEN1 (7H8) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human FEN-1 protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


WB, ICC/IF

Dilutions


WB: 1/1000

;ICC: 1/400

Validations



Western blot detection of FEN-1 in HeLa, Jurkat, 3T3, COS7, PC-12, C6, Raw264.7, Raji and CHO-K1 cell lysates using FEN-1 mouse mAb (1:1000 diluted). Predicted band size: 45KDa. Observed band size: 45KDa.



Immunocytochemistry staining of HeLa cells fixed with 4% Paraformaldehyde and using FEN-1 mouse mAb (dilution 1:400).