

PBX4 antibody



Catalog Number: 113616

Product name

PBX4 antibody

Clonality

Polyclonal

Immunogen

[Human PBX4 Recombinant protein \(GST tag & His tag\)](#)

Ig Type

Rabbit IgG

Specificity

Human; other species not tested.

Applications

ELISA, WB

Antibody description

PBX4 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, A2780 cells, COLO 320 cells, Jurkat cells, SKOV-3 cells. Observed molecular weight by Western-blot: 70 kDa

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of PBX4 recombinant protein (Accession Number: NM_025245). Purification method: Antigen Affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113616(PBX4 Antibody) at dilution of 1:600