

Anti-Bcl-2 Rabbit antibody



Catalog Number: 167316

Product name

Anti-Bcl-2 Rabbit antibody

Specificity

Human, Mouse

Antibody description

Bcl-2 Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human Bcl-2

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC

Dilutions


WB: 1/1000

IHC: 1/20-1/500


ICC/IF: 1/20

FC: 1/50

Validations




Western blot detection of Bcl-2 in MCF-7 cell lysates using Bcl-2 Rabbit pAb(1:1000 diluted). Predicted band size: 26 kDa. Observed



band size: 26 kDa.

Western blot detection of Bcl-2 in K562, CHO-K1, 3T3, HeLa cell lysates using Bcl-2 Rabbit pAb(1:1000 diluted). Predicted band size: 26 kDa. Observed band size: 26 kDa.



Immunohistochemistry of Bcl-2 in paraffin-embedded Human breast cancer tissue using Bcl-2 Rabbit pAb at dilution 1/100